

A faster recovery is at your fingertips.




PLEURAFLOW
ACTIVE CLEARANCE TECHNOLOGY

INTRODUCING THE NEXT-GEN PLEURAFLOW.

OUR MOST VERSATILE ACTIVE CLEARANCE
TECHNOLOGY EVER.

510(k)
cleared

The Only Clinically Proven Active Clearance System






PleuraFlow ACT® is more user-friendly than ever. Our next-gen 510K-cleared technology has been optimized to deliver additional power with the push of a button. Giving you the control needed to prevent or minimize the risks of retained blood complications like tamponade and pleural or pericardial effusions.

510(k)
cleared

Ordering Information



Available in a range of sizes and configurations:

TYPE		CODE	SIZE	# OF EYELETS	EFFECTIVE DRAINAGE LENGTH (EDL)
STRAIGHT		PFFG3-20	20FR	6	4" (10.2 cm)
		PFFG3-24	24FR	6	4" (10.2 cm)
		PFFG3-28	28FR	6	4" (10.2 cm)
		PFFG3-32	32FR	6	4" (10.2 cm)
LONG (XDL)		PFFG3-20XDL	20FR	15	10" (25.4 cm)
		PFFG3-24XDL	24FR	15	10" (25.4 cm)
RIGHT ANGLE (RA)		PFFG3-24RA	24FR	6	4" (10.2 cm)
		PFFG3-28RA	28FR	6	4" (10.2 cm)
		PFFG3-32RA	32FR	6	4" (10.2 cm)

Discover how PleuraFlow helps make recovery beautiful:

CUSTOMER SERVICE: 1-714-916-5007

FAX: 1-714-916-5019

TOLL FREE: 1-844-257-3569 (1-844-CLR-FLOW)

E-MAIL: customerservice@clearflow.com



Manufactured For:
ClearFlow, Inc.
16 Technology Drive, Suite 150
Irvine, CA, US, 92618
www.clearflow.com

Copyright © 2021 ClearFlow, Inc. All Rights Reserved.
™ Trademarks and ® Registered Trademarks of
ClearFlow, Inc.

ML296-B